

Iodine Purification

<u>Iodine temp</u>	1 Quart Water Bottle		2.5 Gallon Water Cube	
	<u>Amount</u>	<u>Wait</u>	<u>Amount</u>	<u>Wait</u>
Body temp	5 cc	40 min	2.5 oz	40 min
Backpack	8 cc	40 min	4 oz	40 min
Freezing	10 cc	40 min	5 oz	40 min

Double the amount of iodine to cut the time in half.
Quadruple the amount of iodine in very cold or questionable water.